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Thank you,

EN

for purchasing this S.C.A.T. Europe product. Please read this operating manual carefully before use and especially observe the safety notices! For service and support, please contact our customer service or your local distributor. For a distributor in your country, please contact us on www.scats-europe.com.

Safety Funnels

Thank you for purchasing our (**electrostatic conductive**) **Safety Funnel**. This Funnel has been designed in such a way to comply with the provisions of CENELEC 50404 (Code of practice for the avoidance of hazards due to static electricity) and the German guideline TRGS 727 (formerly TRBS2153) TRGS 727 can be consulted and downloaded on the website www.baua.de. Where to obtain CLC/TR50404 in your country, please request at info@cencenelec.eu

Important notice:

In order to ensure safe operation we strongly recommend to read this operating manual through carefully before putting your Funnel into service and to comply with the advices when using it. Please check the chemical resistance of the Funnel against the used substances. Also observe the regulations concerning environmental protection and safety in force locally.

Content

Overview	18
Funnels with Ball Valve.....	19
Overview	20
Funnels with Hinged Lid.....	21
Overview	22
Funnels with Level Control.....	23
Intended use.....	24
Please notice before any use!	24
Grounding and Antistatic.....	25
Thermal Properties.....	27
Personal Protective Equipment	27
Maintenance and Care	28
Guarantee and Liability	29

**Scan QR Code =
PDF Download**

Overview

EN

Funnels with Ball Valve

The black Funnels consist of electrostatic conductive plastic (polyethylene high-density PE-HD) and are produced due to quality management standard ISO 9001. The Funnels exhibit the following electrostatic properties: The Funnels are electrostatic conductive.

Funnel sizes

(Thread size of the screw cap)

	Part No.	
	Material PE-HD electrostatic conductive (black)	Material PE-HD (white)
C GL 45	117 622	117 642
C S 50	117 629	117 649
C S 51	117 624	117 644
C S 55	117 625	117 645
C S 60/61	117 621	117 641
C S 65	117 626	117 646
C S 70/71	117 628	117 648
C B 83	117 627	117 647
C S 90	117 623	117 643

Replacement parts:

A Sieve	117 620	117 640
B Splash guard	117 631	117 639
D Grounding cable	117 982	-

Overview

EN

B Sieve

A Hinged Lid

Funnel body

C Screw cap

D Lance

E Grounding cable

XL Funnel

20

www.s...europe.com

Funnels with Hinged Lid

The black Funnels consist of electrostatic conductive plastic (polyethylene high-density PE-HD) and are produced due to quality management standard ISO 9001. The Funnels exhibit the following electrostatic properties: The Funnels are electrostatic conductive.

Funnel sizes

Thread size of the screw cap

	Part No.			
	Material PE-HD electrostatic conductive (black)		Material PE-HD (white)	
	with Lance D	without Lance D	with Lance D	without Lance D
C GL 45	-	118 962	118 992	118 952
C GL 45 XL Funnel	-	117 633	-	-
C S 50	118 985	-	118 995	118 955
C S 51	118 983	-	118 993	118 953
C S 55	118 981	118 961	118 991	118 951
C S 60/61	118 980	118 960	118 990	118 950
C S 60/61 XL Funnel	-	117 634	-	-
C S 65	118 984	118 964	118 994	118 954

Replacement parts:

A Hinged Lid	118 988	118 998
B Sieve	118 989	118 999
E Grounding cable	117 982	-

Overview

EN

Funnels with Level Control

The black Funnels consist of electrostatic conductive plastic (polyethylene high-density PE-HD) and are produced due to quality management standard ISO 9001. The Funnels exhibit the following electrostatic properties: The Funnels are electrostatic conductive.

	<p>Part No.</p>						
<p>Thread size of the screw cap</p> <p>C S 60/61</p>	<p>Material PE-HD electrostatic conductive (black)</p> <p>with Level Control</p> <p>117 637</p>						
	<p>Replacement parts:</p>						
	<table border="1"> <tr> <td data-bbox="265 896 609 950">A Hinged Lid</td> <td data-bbox="609 896 883 950">118 988</td> </tr> <tr> <td data-bbox="265 950 609 998">B Sieve</td> <td data-bbox="609 950 883 998">on request</td> </tr> <tr> <td data-bbox="265 998 609 1049">E Grounding cable</td> <td data-bbox="609 998 883 1049">117 982</td> </tr> </table>	A Hinged Lid	118 988	B Sieve	on request	E Grounding cable	117 982
A Hinged Lid	118 988						
B Sieve	on request						
E Grounding cable	117 982						
<p>Accessories:</p>	<table border="1"> <tr> <td data-bbox="265 1144 609 1192">F S 60/61 Angle Adapter</td> <td data-bbox="609 1144 883 1192">108 058</td> </tr> </table>	F S 60/61 Angle Adapter	108 058				
F S 60/61 Angle Adapter	108 058						

Place laboratory glass and sample bottles on the sieve and let them drip of easily and safely. Also made of electrostatic conductive PE-HD.

Intended use

S.C.A.T. Europe safety funnels are destined for collecting fluid chemicals, such as solvents, solvent mixtures, water and suspensions. The chemicals may be electrostatic conductive, anti-static or insulating. Any usage beyond this intended use for safety reasons is not allowed.

Funnels with Ball Valve: By disposing waste with small particles please rinse enough to ensure the proper function of the ball valve.

Please notice before any use!

When working with the safety funnel, please wear appropriate protective clothing! Always follow the local safety regulations of your company and ask your safety officer in case of any questions!

Please make sure that the safety funnel is installed tightly and free of congestions. When using a safety funnel with ball valve, the ball has to be movable and smooth-running. If the mechanical parts are blocked by dried or viscous fluids, please clean the safety funnel well before using (see „Maintenance and Care“).

Please check the chemical resistance of the funnel material (PE-HD) against the used chemicals. Due to the variety and different constitution of the chemicals available, we do not assume any guarantee or liability regarding chemical resistance. The applicability has to be checked by the operator in each case. More information on chemical resistance can be provided by the manufacturer of the respective chemicals and can be found in the download area of www.scats-europe.com. Please regard the safety instructions and regulations given there.

Use the safety funnel only in case of very good chemical resistance. „Very good“ means that the PE-HD material may not show any damage after 30 days of exposure to the chemical.

Do not use the safety funnel in case of insufficient chemical resistance of the PE-HD material. This may cause damage and loss of the electrostatic conductive characteristics.

Always follow the local safety regulations of your company and ask your safety officer in case of any questions!

Grounding and Antistatic

When working with highly flammable fluids with low flashpoint, any device made of electrostatic conductive materials must be grounded safely. Any electrostatic conductive safety funnels made of electrostatic conductive PE-HD (black) must be grounded actively during operation.

For grounding, use only the provided grounding cable.

Warning: The installation has to be executed by specialized people only! Please notice the internal safety advises of your company.

The safety funnel is deemed to be grounded if the resistance to ground is max. 10^6 Ohm at a voltage of 500V. If the grounding has been carried out correctly, no electrostatic charges are due when handling electroconductive or antistatic fluids.

When operating with suspensions and insulating fluids, an accumulation of electrostatic charge is possible. This accumulation develops autonomously, independent from the safety funnel. Accumulation on the safety funnel itself may not occur if the grounding has been carried out correctly.

The following factors can affect the electrostatic characteristics of electrostatic conductive safety funnels negatively:

- **The fluid or suspension is sticky and/or insulating.** This may create a insulating coating at the inner surface of the safety funnel.
- **Diffusion of an insulating fluid** (typically non-conductive solvents) through the plastic material. Over a longer period of time, this may cause loss of the safety funnel's electrostatic conductive characteristics.

Important notice: Electrostatic characteristics of the PE-HD material may change over a longer period of time. The resistance to the ground has to be checked at least once a year. The resistance to ground has to be max. 10^6 Ohm at a voltage of 500V.

In case of deterioration, the safety funnel may no longer be regarded as antistatic and has to be replaced.

Please always follow the local safety regulations of your company and ask your safety officer in case of any questions!

Thermal Properties

The permanent operating temperature of PE without mechanical strain is $-50\text{ }^{\circ}\text{C}$ to $80\text{ }^{\circ}\text{C}$. To avoid the funnel from deformations you should not fill in liquids with a temperature higher than $60\text{ }^{\circ}\text{C}$.

EN

Personal Protective Equipment

Always wear appropriate protective clothing to protect yourself against spillings and vapours during operation.

During maintenance and cleaning of the safety funnel, always wear protective gloves and glasses!

Always follow the local safety regulations of your company and ask your safety officer in case of any questions!

Maintenance and Care

EN

When working with the safety funnel and chemicals, always wear protective clothing and follow the local safety regulations of your company! Ask your safety officer in case of any questions!

After each use, the safety funnel has to be cleaned and purged with an appropriate solvent to avoid congestions or blocking of the ball valve. In-depth cleaning avoids non-conductive coatings or diffusion that may affect the electrostatic conductive characteristics of the material surface and the anti-static performance of the safety funnel.

Please check the movability of the ball valve regularly. Always purge the funnel after using sticky or viscous fluids. Remove the sieve and the splash guard, and clean them separately.

Purge funnels regularly, e.g. when changing containers, or earlier where necessary. Suitable solvents for purging are water, acetonitrile or isopropyl alcohol. Clean especially movable parts regularly and check them for error-free operation, e.g. the floater when using funnels with level control. This makes sure that there will be no jamming and crystallization caused by drying due to the variety of possible waste liquids.

Please always use wet cloth for external cleaning. Dry rubbing may cause static charge and sparking when components are not grounded sufficiently.

Guarantee and Liability


Please do not make any technical changes or modifications. If the safety funnel is technically modified by the operator, S.C.A.T. Europe GmbH does not assume any guarantee or liability.

If you need a technically modified product for a special application, please contact us. Our product specialists will help you in developing a customized solution.

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User Manual

PDF Download

The latest version of **User Manual** relating to your product is always available for you to download. Please visit www.scat-europe.com or scan this **QR-Code** using the camera of your Smartphone/Tablet.

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A stylized illustration of a hand with a pink and white polka-dot pattern, wearing a blue sleeve, holding a white rectangular sign. The sign contains text about finding a local distributor. A QR code is located at the bottom right of the sign.

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© S 60/61