



**WTW Safety Data Sheet**  
According to EC Directives 91/155/EEC

**1. Identification of the substance/preparation and of the company/undertaking**

*Product name:*

**Standard B      Meas. Range 1      1,6 mg/l PO4-P**

**Ord-No.:            821 121**  
**Model:             S/BI-P 211**

*Manufacturer:*

WTW Wissenschaftlich-Technische Werkstätten GmbH, Dr.-Karl-Slevogt-Straße 1, D-82362 Weilheim, Germany, Tel.: (0881) 183-0, Fax.: (0881) 62539

Emergency telephone:  
(001)-703-527-3887

**2. Composition/information on ingredients**

*Chemical characterization:*            Aqueous saline solution

*Ingredients:*                                Potassium dihydrogen phosphate

*Hazard symbols:*                        No hazardous product  
*R-phrases:*                                /  
*Content:*                                    < 0,001 %  
*CAS-No.:*                                  7778-77-0  
*EINECS-No.:*                              231-913-4

**3. Hazards identification**

No hazardous product.

**4. First aid measures**

*After inhalation:*                        Fresh air  
*After skin contact:*                    Wash off with plenty of water.  
Remove contaminated clothing.  
*After eye contact:*                    Rinse out with plenty of water  
with the eyelid held wide open.  
*After ingestion:*                        Drink plenty of water.  
Consult doctor if feeling unwell.

## WTW Safety Data Sheet

### Identification of the substance/preparation

Standard B      Meas. Range 1      1,6 mg/l PO4-P

---

#### 5. Fire-fighting measures

*Suitable extinguishing media:*      To suit environment  
*Unsuitable extinguishing media:*      No information available.  
*Special risks:*      None  
*Other information:*      Non combustible

#### 6. Accidental release measures

*Personal-related precautionary measures:*      None  
*Environmental-protection measures:*      None  
*Procedures for cleaning/absorption:*      Clean up affected area with water.

#### 7. Handling and storage

*Handling:*      No further requirements  
*Storage:*      Store tightly closed.

#### 8. Exposure control/personal protection

*Specific control parameters:*      None  
Personal protective equipment  
*Respiratory protection:*      Not required  
*Eye protection:*      Protective glasses  
*Hand protection:*      Gloves  
*Industrial hygiene:*      Wash hands after working with substance.



# WTW Safety Data Sheet

## Identification of the substance/preparation

Standard B      Meas. Range 1      1,6 mg/l PO4-P

---

### 13. Disposal considerations

*Products:*                      There are no uniform regulations for the disposal of chemicals or residues. The disposal is regulated through corresponding regional laws and regulations. We recommend that you contact either the authorities in charge or approved waste disposal companies which will advise you on how to dispose the special waste.

*Packaging:*                    Handle contaminated packaging in the same way as the substance itself. If not officially specified differently, noncontaminated packaging may be treated like household waste or recycled.

### 14. Transport information

Not subject to transport regulations.

### 15. Regulatory information

Labelling according to EC directives

*Symbol:*                                      No hazardous product  
*R-phrases:* /  
*S-phrases:* /

German regulations

*Water pollution class;*                      Nonpolluting substance

*Regulation on combustible liquids;* VbF: Non combustible

*Other regulations:*                              None

### 16. Other information

Date of issue:                                  18.2.2004  
Supersedes edition of                      20.8.2001  
Reason for alteration:                        Generall update

The information contained herein is based on the present state of our knowledge. It characterizes the product with regard to the appropriate safety precautions. It does not represent a guarantee of the properties of the product.