



E 404/1








Mobile unit
for the optimum loading of up to 38 pipettes.

- Suitable for 38 pipettes in 3 rows
- 1st Row: 10 pipettes x 100 ml (up to 550 mm), bar distance 20 mm
- 2nd Row: 14 pipettes x 25 ml, bar distance 26 mm
- 3rd Row: 14 pipettes x 10 ml, bar distance 26 mm
- Without bracket for magnetic strip for mobile unit recognition



Product details - E 404/1

EAN: 4002514087061 / Article no. 69540402D / Mat. no. 06727910

Product category 	
Washer-disinfectors, laboratory	•
Machine/accessory categorisation 	
PG 8535	•
PG 8536	•
Product category 	
Mobile unit for pipettes	•
Product properties 	
Material	Stainless steel Plastic
Colour	Stainless steel Black
Number of rows	3
Area of application 	
Reprocessing pipettes/vials	•
Capacity 	
Pipettes [number]	38
Dimensions and weight 	
External dimensions, height in mm	502
External dimensions, width in mm	535
External dimensions, depth in mm	517
Net weight in kg	5,700

Subject to technical changes; no liability accepted for the accuracy of the information given!

Distributed by:



Tallaght Business Park
Whitestown, Dublin 24,
Ireland
D24 RFK3

Tel: (01) 4523432
Fax: (01) 4523967
E-mail: info@labunlimited.com
Web: www.labunlimited.com

Quatro House, Frimley Road,
Camberley,
United Kingdom
GU16 7ER

Tel: 08452 30 40 30
Fax: 08452 30 50 30
E-mail: info@labunlimited.co.uk
Web: www.labunlimited.co.uk