


# SoftGrip™ Pipette Replacement Parts

SoftGrip™ pipettes can be easily serviced and calibrated in the lab. If you prefer to service your own pipette, replacement parts are readily available to ensure long pipette life. Service and calibration are also available through Hamilton; see page 12 for information on Hamilton's Complete Pipette Service.


## 2 µL Pipette



### SoftGrip Single Channel Pipette, 2 µL

Part Description	2 µL Part No.
Plunger Cap	54190-02
Lower Body Assembly**	4118-07
Seal Spring	1288-12
Seal Spacer	4029-01
Seal Kit*	1708-79
2 µL Tip	4028-01
Ejector Sleeve	54758-01

## 5 µL - 300 µL Pipettes



### SoftGrip Single Channel Pipette, 5 to 300 µL


Part Description	Part Numbers							
	5 µL	10 µL	25 µL	50 µL	100 µL	200 µL	250 µL	300 µL
Plunger Cap	54190-13	54190-11	54190-03	54190-06	54190-12	54190-07	54190-08	54190-10
Spring Spacer	34189 for 5-300 µL pipettes							
Seal Spring	1390-01A for 5-300 µL pipettes							
Seal Spacer	34106	34106	34136	34062	34107	34012	34012	34012
Seal Tube	34064	34064						
Seal Kit*	1708-70	1708-70	1708-71	1708-72	1708-73	1708-74	1708-75	1708-76
Lower Body Assembly**	4118-08	4118-09	4118-10	4118-11	4118-12	4118-13	4118-14	4118-15
Ejector Sleeve	54758-01	54758-01	54757-01	54757-01	54757-01	54757-01	54757-01	54757-01

\* Seal Kit contains 5 PTFE seals with O-rings, and 1 tube of lubricant.

\*\* Lower body assembly comes with one seal kit.

# SoftGrip™ Pipette Replacement Parts

## 500 µL - 1 mL Pipettes




### SoftGrip Single Channel Pipette, 500 to 1000 µL

Part Description	Part Numbers	
	500 µL	1000 µL
Plunger Cap	54190-04	54190-01
Retainer Assembly	34196	34196
Eject Spring	1288-06	1288-06
Adapter	34154	34154
Seal Spring	1390-01A	1390-01A

Part Description	Part Numbers	
	500 µL	1000 µL
Seal Spacer	34132	34133
Seal Kit*	1708-77	1708-78
Lower Body Assembly**	4118-16	4118-17
Ejector Sleeve	54756-01	54756-01

## 50 µL and 300 µL Multi-Channel Pipettes



### SoftGrip Multi-Channel Pipette

Part Description	Part Numbers			
	8 x 50 µL	12 x 50 µL	8 x 300 µL	12 x 300 µL
Plunger Cap	54190-06	54190-06	54190-10	54190-10
Spring Clip	34354	34354	34354	34354
Piston Chamber Assembly	6385-01	6385-01	6385-02	6385-02
Ejector Plate	4923-01	4922-01	4923-01	4922-01
Cover	4917-01	4916-01	4917-01	4916-01

\* Seal Kit contains 5 PTFE seals with O-rings, and 1 tube of lubricant.

\*\* Lower body assembly comes with one seal kit.

Distributed by:



**Carl Stuart Limited**

ADVANCED APPLIED TECHNOLOGIES

Contact Us:

Irl Ph: 01 4523432

UK Ph: 08452 30 40 30

Web: [www.carlstuart.com](http://www.carlstuart.com)

Email: [info@carlstuart.com](mailto:info@carlstuart.com)

SoftGrip™  
Pipettes

SoftTouch™  
Electronic Pipettes

AdvantiTip™  
Pipette Tips

SoftAide™  
Pipette Controller

SoftTop™  
Quik Dispenser

Replacement Parts  
and Accessories