

# Jeio Tech 198 Litre Stackable Incubator Shaker (ISS-7200)

**Product Code:** AAH239132K

**Brand:** Jeio Tech

Maximise your laboratory space with the Jeio Tech 198 Litre Stackable Incubator Shaker (ISS-7200).

This large, multi-purpose incubator shakers allows the user to freely organise different sample sizes and can stack up to three units. It can accommodate, for example, six 4 litre flasks and also thirty-nine 250ml flasks. If you decide to stack all three units, you can triple the capacity, although the highest rpm will be reduced.

The incubator shaker's door is opened upwards and stays open for hands-free handling. A slide type universal platform will let you easily access all of your samples, even ones placed at the back.

An intuitive touch-screen LCD controller with easy icons and logical menus is included. Universal platform is required (not included).

## Features

### Performance

- **Temperature Range:** Ambient +10°C to 80°C
- Competitive high shaking speed up to 400rpm (stackable up to 200rpm).
- Microprocessor PID control
- Three point temperature calibration/automatic tuning
- Wide speed range even with heavy workload. Ideal DD (Direct Drive) brushless motor and triple cam system support heavy load, quiet, maintenance-free operation and deliver high torque at high speeds and smooth rotation
- Best effort run function intelligently manages its rpm to keep shaking even workload is out of its capacity. Automatic shaking speed adjustment in case of excessive workload such as unbalanced load placement, unusual vibrations caused from unstable floor or external shock
- Pleasant test environment by the smooth acceleration and deceleration control. Smooth start and stop function prevents chemical spills from flasks or test tubes
- High-velocity fan ensures uniform temperature distribution and rapid thermal recovery after the door is opened
- Air-tight silicone door seal to provide excellent temperature uniformity

### Convenience

- Intuitive touch-screen LCD controller with easy icons and logical menus
- Easy-set digital timer for shaking operation (1 min. to 999 hr. 59 min.)
- Wide sample monitoring. Clear observation of samples without affecting inner chamber's environment thanks to the transparent viewing window with bright LED lamps
- Repetitive tasks can be easily performed as the platform stops where it starts. Stopping the shaking platform always at the same position is highly beneficial for automated closing or sampling process
- Built-in electrical outlet with a safety cover inside the chamber
- Corrosion resistance stainless steel interior
- Easy-access drain system with a quick disconnect valve offers great convenience to clean up spills
- Built-in RS-232 and USB port for external control and data collection
- The unit comes with retractable foot caster, beneficial for easy mobility during installation or relocation

### Safety

- Automatic run after power interruption
- Temperature and shaking speed deviation alarm
- Over-current protection
- **Triple Independent Temperature Monitor System:** Temperature deviation alarm (high/low), electronic temperature limiter (high/low), and mechanical over-temperature limiter (high)



# Specifications

<b>Volume (L)</b>			198
<b>Temperature Range (°C)</b>			Amb. +10 to 80
<b>Fluctuation at 37°C</b>			± 0.1
<b>Variation at 37°C</b>			± 0.5
<b>Orbit Diameter (mm)</b>			50.8
<b>Speed Range (rpm)</b>			30 to 300
<b>Speed Range, stacked (rpm)</b>			30 to 200
<b>Timer</b>			1 min. to 999 hr. 59 min.
<b>Max. Load (kg)</b>			20 at 300rpm
<b>Max. Load, stacked (kg)</b>			15 at 200rpm 20 at 150rpm
<b>Platform Dimensions (W x D)</b>			755 x 520mm
<b>Internal Dimensions (W x D x H)</b>			809 x 624 x 393mm
<b>External Dimensions (W x D x H)</b>			1310 x 996 x 755mm
<b>Weight (kg)</b>			365
<b>Max. Mountable Flask Clamps</b>	<b>Stainless Steel</b>	<b>For 50ml</b>	88
		<b>For 100ml</b>	68
		<b>For 250ml</b>	39
		<b>For 300ml</b>	
		<b>For 500ml</b>	28
		<b>For 1 Litre</b>	18
		<b>For 2 Litres</b>	11
		<b>For 2.8 Litres</b>	8
		<b>For 4 Litres</b>	6
	<b>Plastic</b>	<b>For 50ml</b>	88
		<b>For 100 - 125ml</b>	48
		<b>For 200ml</b>	33
		<b>For 250ml</b>	
		<b>For 300ml</b>	
		<b>For 500ml</b>	24
		<b>For 1 Litre</b>	12
<b>For 2 Litres</b>	8		

**Lab Unlimited**  
CARL STUART GROUP

A CALIBRE SCIENTIFIC COMPANY

#### Carl Stuart Limited

Tallaght Business Park  
Whitestown  
Dublin 24  
D24 RFK3  
Ireland

Ph: 00353 1 4523432  
Fax: 00353 1 4523967  
E-mail: [onlinestore@labunlimited.com](mailto:onlinestore@labunlimited.com)

#### Carl Stuart UK Limited

Calibre Scientific UK  
Unit 5A, R-Evolution @  
The Advanced Manufacturing Park  
Selden Way  
Rotherham  
S60 5XA  
United Kingdom

Ph: 08452 30 40 30  
Fax: 08452 30 50 30  
E-mail: [onlinestore@labunlimited.co.uk](mailto:onlinestore@labunlimited.co.uk)